

+ ALLEVYN[◇] GENTLE Foam Dressing is a silicone non-bordered dressing designed for when you need some gentle adhesion to give you a helping hand as you apply dressings on delicate areas of the body.^{1,2}

Smith+Nephew

ALLEVYN[◇] GENTLE
Foam Dressing

Application



- 1**
- Remove first liner and anchor the adhesive side of the dressing to the skin



- 2**
- Remove remaining liner/liners



- 3**
- Smooth over
 - Ensure dressing covers the wound
 - Do not stretch the dressing
 - The dressing can be repositioned as required



- 4**
- Secure with secondary retention e.g. OPSITE[◇] FLEXIFIX Transparent Film Roll, OPSITE FLEXIFIX GENTLE Transparent Silicone Adhesive Film Roll, tape, bandage or compression therapy


Removal

Lift one corner and slowly peel back until completely removed.

Frequency of change

During the early stages of treatment inspect the dressing frequently. Dressings can be left in place for up to 7 days. Dressings should be changed depending on the condition of the wound and surrounding skin or when exudate is within 1.5cm (5/16in) of the dressing edge.

The diagram below is for guidance only and the decision of when to change should be dependent upon clinical assessment and local protocols under the guidance of a Healthcare Professional.



ALLEVYN GENTLE Foam Dressing ordering codes

Size	Code	Carton	Size	Code	Carton
5cm x 5cm / 2in x 2in	66802128	10	15cm x 15cm / 6in x 6in	66802131	10
10cm x 10cm / 4in x 4in	66802129	10	20cm x 20cm / 8in x 8in	66802132	10
10cm x 20cm / 4in x 8in	66802130	10	20cm x 50cm / 8in x 20in	66802133	2

Advanced Wound Management
Smith & Nephew Medical Ltd
Croxley Park, Building 5,
Hatters Lane, Watford,
Hertfordshire WD18 8YE, UK

www.smith-nephew.com
T +44 (0) 1923 477100
F +44 (0) 1923 477101

◇Trademark of Smith+Nephew
All Trademarks acknowledged
© February 2021 Smith+Nephew
AWM-AWC-29139-uki

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

References

1. Smith+Nephew 2014. An open, prospective, randomised volunteer trial to assess the retention qualities of 4 non-bordered silicone foam dressings under secondary retention. Internal report. VTR/HVT079/ST965. 2. Vowden K. An open, prospective randomised, multi-centre clinical evaluation of a hydrocellular silicone foam dressing* in the management of exuding chronic and acute wounds. Paper presented at: EWMA; 2011; Brussels.