

Medi UK introduces the world's fastest ABPI

November 2017

medi know a thing or two about effective compression – we also know that regular ABPI checks are required to ensure that compression is safe to apply. Which is why we are delighted to bring you the quickest and simplest method of checking your patients' ABPI – the **ABPI MD**

- Achieve CQUIN targets easily by taking accurate ABPI measurements in one minute.

Early diagnosis of Peripheral Arterial Disease (PAD) in primary healthcare is crucial. Over 70% of patients do not know about the disease because they do not feel or recognise the symptoms until severe complications occur. Ankle Brachial Pressure Index (ABPI) measurement is a simple solution for diagnosing PAD. It compares blood pressure in the arm and legs. With the new ABPI MD device, the procedure is reliable, objective and even possible to be performed as a screening device in primary healthcare.

Features:

- No interpretation of results required – results displayed show ABPI, BP and Pulse readings
- No need to rest patient beforehand
- ABPI results and brachial blood pressure in just 1 minute
- Simultaneous measurement of left and right ABPI
- Unique algorithm for ankle blood pressure calculation
- Simultaneous measurement to avoid error from blood pressure drift
- Cuff-based technology sensors detect smallest changes in volume
- Advanced error detection system prevents false reporting
- Data can be saved and detailed reports can be printed

Revolutionise your clinic now!

For more information, email admin.team@mediuk.co.uk or Tel. 01432 373500.

<https://www.mediuk.co.uk/products/mesi-abpi-md/>

Watch the short video here - <https://player.vimeo.com/video/237564652>

medi UK Ltd, Plough Lane, Hereford, HR4 0EL.